

School Code: 10908



CBSE Affiliation No. : 1730578

BHAGAT PUBLIC SR. SEC. SCHOOL

ALANIYA, KOTA



CHEMISTRY PRACTICAL

MAJOR EXPERIMENT

2020-21

CLASS - XII



ADDRESS

CAMPUS : N.H. 12, JHALAWAR ROAD, ALANIYA, KOTA-325003, PH :0744-2832113, 9649991123


Email : bpssschool@gmail.com

2


Experiment No: 1

AIM:- To prepare 250 ml of M/20 standard solution of F.A.S. Using this solution find out the molarity and strength of the given KMnO_4 solution.

Theory :- KMnO_4 oxidises Fe^{2+} ions into Fe^{3+} ions in acidic medium in cold and itself reduced to colourless Mn^{2+} ions



Or



Molecular Mass of $\text{KMnO}_4 = 158$

Molecular Mass of F.A.S. = 392

Mass of required FAS for preparing 250 ml of

$$\text{M/20 solution} = \frac{392}{20} \times \frac{250}{1000} = 4.9 \text{ gm}$$

Apparatus:- Chemical balance, weight box, watch glass, measuring flask, funnel, burette, pipette, conical flask, tile, burette stand etc.

Materials- F.A.S Crystals (Mohr's Salt), Dilute H_2SO_4 , Distilled H_2O , KMnO_4 solution.

Indicator - KMnO_4 solution act as a self indicator.

End point - colourless to Permanent pink

Observation -

- (i) Mass of watch glass (a) = 18.6192 gm
- (ii) Mass of mohr's salt required (b) = 4.9 gm
- (iii) Mass of (a) + (b) = 23.5192 gm
- (iv) Volume of mohr's salt solution prepared = 250 ml
- (v) Molarity of mohr's salt solution = M/20

Observation table -

S.N.	Volume of F.A.S. Solution (V_1 ml)	Reading of KMnO_4 solution			Concordant volume of KMnO_4 (V_2 ml)
		Initial (a) (ml)	final (b) (ml)	Used (b-a) (ml)	
1	20 ml	0.0	19.8	19.8 ml	
2	20 ml	0.0	19.6	19.6 ml	19.6 ml
3	20 ml	0.0	19.6	19.6 ml	

Calculations- (A) Molarity of KMnO_4 solution (M_2) = ? .
Using formula $M_1V_1 = 5 M_2V_2$

$$M_2 = \frac{M_1V_1}{5 \times V_2} = \frac{1}{20} \times \frac{20}{5 \times 19.6}$$

$$M_2 = 0.0102 \text{ M}$$

$$M_1 = \text{M/20}$$

$$V_1 = 20 \text{ ml}$$

$$M_2 = ?$$

$$V_2 = 19.6 \text{ ml}$$

(due to 2 moles of KMnO_4 reacts with 10 moles of FAS in overall balance equation)

(B) Strength of KMnO_4 Solution in gm /lit = ?
= molarity (M_2) \times molecular mass
= 0.0102×158
= 1.6116 gm /lit

Result – (1) Molarity of the given KMnO_4 solution = 0.0102 M
(2) Strength of the given KMnO_4 solution = 1.6116 gm/lit .


Precaution

- (i) All the volumetric apparatus should be washed well before use.
- (ii) Rinse burette with the solution of KMnO_4 and pipette with the solution of FAS. And wash the titration flask with distilled H_2O after every titration.
- (iii) Always read the upper meniscus in case of coloured solutions.
- (iv) Always use freshly prepared KMnO_4 as it decomposed on keeping.
- (v) Excess of dilute H_2SO_4 should be add in FAS.


Experiment No : 2

AIM:- To prepare 250 ml of M/30 standard solution of F.A.S. Using this solution find out the molarity and strength of the given KMnO_4 solution.

Theory :- KMnO_4 oxidises Fe^{2+} ions into Fe^{3+} ions in acidic medium in cold and itself reduced to colourless Mn^{2+} ions



Or



Molecular Mass of $\text{KMnO}_4 = 158$

Molecular Mass of F.A.S. = 392

Mass of required FAS for preparing 250 ml of

$$\text{M/30 solution} = \frac{392}{30} \times \frac{250}{1000} = 3.26 \text{ gm}$$

Apparatus:- Chemical balance, weight box, watch glass, measuring flask, funnel, burette, pipette, conical flask, tile, burette stand etc.

Materials:- F.A.S Crystals (Mohr's Salt), Dilute H_2SO_4 , Distilled H_2O , KMnO_4 solution.

Indicator – KMnO_4 solution act as a self indicator.

End point – colourless to Permanent pink

- Observation –**
- (i) Mass of watch glass (a) = 18.5974 gm
 - (ii) Mass of mohr's salt required (b) = 3.26 gm
 - (iii) Mass of (a) + (b) = 21.8574 gm
 - (iv) Volume of mohr's salt solution prepared = 250 ml
 - (v) Molarity of mohr's salt solution = M/30

Observation table –

S.N.	Volume of F.A.S. Solution (V_1 ml)	Reading of KMnO_4 solution			Concordant volume of KMnO_4 (V_2 ml)
		Initial (a) (ml)	final (b) (ml)	Used (b-a) (ml)	
1	20 ml	0.0	19.9	19.9 ml	
2	20 ml	0.0	19.7	19.7 ml	19.7 ml
3	20 ml	0.0	19.7	19.7 ml	

Calculations- (A) Molarity of KMnO_4 solution (M_2) = ?

Using formula

$$M_1V_1 = 5 M_2V_2$$

$$M_2 = \frac{M_1V_1}{5 \times V_2} = \frac{1}{30} \times \frac{20}{5 \times 19.7}$$

$$M_2 = 0.0067 \text{ M}$$

$$M_1 = M/30$$

$$V_1 = 20 \text{ ml}$$

$$M_2 = ?$$

$$V_2 = 19.7 \text{ ml}$$

(due to 2 moles of KMnO_4 reacts with 10 moles of FAS in overall balance equation)

(B) Strength of KMnO_4 Solution in gm /lit = ?

$$= \text{molarity } (M_2) \times \text{molecular mass}$$

$$= 0.0067 \times 158$$

$$= 1.0586 \text{ gm /lit}$$

Result - (1) Molarity of the given KMnO_4 solution = 0.0067 M

(2) Strength of the given KMnO_4 solution = 1.0586 gm/lit.


Precaution

- (i) ii the volumetric apparatus should be washed well before use.
- (ii) Rinse burette with the solution of KMnO_4 and pipette with the solution of FAS. And wash the titration flask with distilled H_2O after every titration.
- (iii) Always read the upper meniscus in case of coloured solutions.
- (iv) Always use freshly prepared KMnO_4 as it decomposed on keeping.
- (v) Excess of dilute H_2SO_4 should be add in FAS.


Experiment No. 3

AIM:-
To prepare 250 ml of M/20 standard solution of Oxalic acid.
Using this solution, find out the Molarity and strength of the given KMnO_4 solution.

Theory – KMnO_4 Oxidises Oxalic acid into CO_2 in acidic medium at a temperature around 60°C and itself get reduced to colourless Mn^{2+} ions



OR



Molecular mass of $\text{KMnO}_4 = 158$

Molecular mass of Oxalic acid = 126

Mass of required Oxalic acid for preparing 250 ml of

$$\text{M/20 solution} = \frac{126}{20} \times \frac{250}{1000} = 1.575 \text{ gm}$$

Apparatus:- Chemical balance, weight box, watch glass, measuring flask, funnel, burette, pipette, conical flask, tile, burette stand etc.

Materials- Hydrated oxalic acid Crystals (Mohr's Salt), dilute H_2SO_4 , distilled H_2O , KMnO_4 solution.

Indicator – KMnO_4 solution act as a self indicator.

End point – Colourless to Permanent pink

Observation –

- Mass of watch glass (a) = 18.6202 gm
- Mass of required oxalic acid (b) = 1.575 gm
- Mass of (a) + (b) = 20.1952 gm
- Volume of oxalic acid solution prepared = 250 ml
- Molarity of oxalic acid solution = M/20

Observation table –

S.N.	Volume of oxalic acid Solution (V_1 ml)	Reading of KMnO_4 solution			Concordant volume of KMnO_4 (V_2 ml)
		Initial (a) (ml)	final (b) (ml)	Used (b-a) (ml)	
1	20 ml	0.0	19.7	19.7 ml	
2	20 ml	0.0	19.5	19.5 ml	19.5 ml
3	20 ml	0.0	19.5	19.5 ml	

Calculations- (A) Molarity of KMnO_4 solution (M_2) = ?

Using formula

$$2M_1V_1 = 5M_2V_2$$

$$M_2 = \frac{2M_1V_1}{5V_2}$$

$$M_2 = \frac{2}{20} \times \frac{20}{5 \times 19.5}$$

$$M_1 = \text{M/20}$$

$$V_1 = 20 \text{ ml}$$

$$M_2 = ?$$

$$V_2 = 19.5 \text{ ml}$$

$$M_2 = 0.0205 \text{ M}$$

(due to 2 mole of KMnO_4 reacts with 5 moles of oxalic acid in overall balance equation)

(B) Strength of KMnO_4 Solution in gm/lit

$$\begin{aligned} &= \text{Molarity } (M_2) \times \text{molecular mass.} \\ &= 0.0205 \times 158 \\ &= 3.2390 \text{ gm/lit} \end{aligned}$$

Result (A) The Molarity of the given KMnO_4 solution = 0.0205 M
(B) The Strength of the given KMnO_4 solution = 3.2390 gm/lit


Precautions

- (i) All the volumetric apparatus should be washed well before use.
- (ii) Rinse burette with the solution of oxalic acid and pipette with the solution of oxalic acid and wash the titration flask with distilled H_2O after every titration
- (iii) Always read the upper meniscus in case of coloured solution.
- (iv) Always use freshly prepared KMnO_4 as it decomposes on keeping.
- (v) Dilute H_2SO_4 should be added in oxalic acid solution.
- (vi) Heat the conical flask should be slowly on a wire gauze unit the moisture appears on the neck of the flask the temperature of the solution should be between $60^\circ - 70^\circ \text{C}$


Experiment No. 4

AIM:- To prepare 250 ml of M/30 standard solution of Oxalic acid.
Using this solution, find out the Molarity and strength of the given KMnO_4 solution.

Theory – KMnO_4 Oxidises Oxalic acid into CO_2 in acidic medium at a temperature around 60°C and itself get reduced to colourless Mn^{2+} ions



OR



Molecular mass of $\text{KMnO}_4 = 158$

Molecular mass of Oxalic acid = 126

Mass of required Oxalic acid for preparing 250 ml of

$$\text{M/30 solution} = \frac{126}{30} \times \frac{250}{1000} = 1.05 \text{ gm}$$

Apparatus:- Chemical balance, weight box, watch glass, measuring flask, funnel, burette, pipette, conical flask, tile, burette stand etc.

Materials- Hydrated oxalic acid Crystals (Mohr's Salt), dilute H_2SO_4 , distilled H_2O , KMnO_4 solution.

Indicator – KMnO_4 solution act as a self indicator.

End point – Colourless to Permanent pink

Observation –

- (i) Mass of watch glass (a) = 18.6192. gm
- (ii) Mass of required oxalic acid (b) = 1.05 gm
- (iii) Mass of (a) + (b) = 19.6692 gm
- (iv) Volume of oxalic acid solution prepared = 250 ml
- (v) Molarity of oxalic acid solution = M/30

Observation table –

S.N.	Volume of oxalic acid Solution (V_1 ml)	Reading of KMnO_4 solution			Concordant volume of KMnO_4 (V_2 ml)
		Initial (a) (ml)	final (b) (ml)	Used (b-a) (ml)	
1	20 ml	0.0	19.9	19.9 ml	
2	20 ml	0.0	19.7	19.7 ml	19.7 ml
3	20 ml	0.0	19.7	19.7 ml	

XII CBSE : Chemistry Practical

Calculations- (A) Molarity of KMnO_4 solution (M_2) = ?

Using formula

$$2M_1V_1 = 5M_2V_2$$

$$M_2 = \frac{2M_1V_1}{5V_2}$$

$$M_2 = \frac{2}{30} \times \frac{20}{5 \times 19.7}$$

$$M_1 = M/30$$

$$V_1 = 20 \text{ ml}$$

$$M_2 = ?$$

$$V_2 = 19.7 \text{ ml}$$

$$M_2 = 0.0135 \text{ M}$$

(due to 2 mole of KMnO_4 reacts with 5 moles of oxalic acid in overall balance equation)

(B) Strength of KMnO_4 Solution in gm/lit

$$= \text{Molarity } (M_2) \times \text{molecular mass.}$$

$$= 0.0135 \times 158$$

$$= 2.1330 \text{ gm/lit}$$

Result (A) The Molarity of the given KMnO_4 solution = 0.0135 M

(B) The Strength of the given KMnO_4 solution = 2.1330 gm/lit

Precautions

- (i) All the volumetric apparatus should be washed well before use.
- (ii) Rinse burette with the solution of oxalic acid and pipette with the solution of oxalic acid and wash the titration flask with distilled H_2O after every titration.
- (iii) Always read the upper meniscus in case of coloured solution.
- (iv) Always use freshly prepared KMnO_4 as it decomposes on keeping.
- (v) Dilute H_2SO_4 should be added in oxalic acid solution.
- (vi) Heat the conical flask should be slowly on a wire gauze unit the moisture appears on the neck of the flask; the temperature of the solution should be between $60^\circ - 70^\circ \text{C}$.

Experiment No. 5

Exercise – Identify the functional group present in the given organic compound.

Physical Properties –


1. Physical State - Liquid
2. Colour – Colour less
3. Odour – Smell like vinegar
4. Solubility in water – Soluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound - Put a crystal of the compound on moist blue litmus paper	Litmus paper turns red	Organic compound is acidic in nature may be -COOH or Ar-OH Functional group.

S.N	Experiment	Observations	Inference
2	Burning Test – Take some organic compound on copper wire and bring it in to the flame.	The compound burns with non-sooty blue flame	Compound is aliphatic.
3	Detection of element (1) Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared)	Green ppt	N-absent
	(ii) Add few drops of dil. H ₂ SO ₄	Colour less solution is obtained	
4	Test for functional group Aqueous solution of compound + a pinch of NaHCO ₃	Brisk effervescence colourless odourless gas evolved	Ar-OH group absent- -COOH group in present.
5	Ester test Organic compound + 2-3 drops H ₂ SO ₄ +0.5 ml C ₂ H ₅ OH. Warm	Fruity smell	-COOH group present

Result – Carboxylic acid [-COOH] group is present in the given organic compound.

Structural Formula



Experiment No. 6

Exercise – Identify the functional group present in the given organic compound.

Physical Properties –


1. Physical State - Liquid
2. Colour – Colour less
3. Odour – Pleasant smell like that of spirit
4. Solubility in water – Soluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound-Put a drop of liquid on moist blue litmus paper and moist red litmus paper respectively	No reaction	Organic comp. is neutral.

S.N	Experiment	Observations	Inference
2	Burning Test – Burning a copper wire dipped in the given liquid into the flame.	The compound burns with non- sooty blue flame	Compound is aliphatic.
3	Detection of element-1. Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared)	Green ppt	
	(ii) Add few drops of dil. H ₂ SO ₄	Solution does not turns green or blue.	N-absent
4	Test for functional group		
	Sodium metal test-Organic liquid + Na piece	Brisk effervescence with the evolution of H ₂ gas	-OH (alcoholic) group is present.
5	<u>Ceric ammonium nitrate test</u> A aqueous solution of the compound + ceric ammonium nitrate solution.	Solution turns red	-OH group is present.
6	Ester Test – Compound + solid CH ₃ COONa + few drops of cone. H ₂ SO ₄ and heat.	Fruity smell	-OH group is present.

Result – The given organic compound has alcoholic (-OH) functional group in it.

Structural formula -



Experiment No.7 ✓

Exercise – Identify the functional group present in the given organic compound.

Physical Properties –


1. Physical State – Solid
2. Colour – Light Pink
3. Odour – Pheolic smell
4. Solubility in water – Soluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound-Put a drop of liquid on moist blue litmus paper	Litmus paper turns red	The given comp. is acidic. May –COOH or Ar-OH group in it.

S.N	Experiment	Observations	Inference
2	Burning Test – Put a crystal of the compound on glass rod and bring it in to the flame	The compound burns with sooty flame	Organic Compound is aromatic.
3	Detection of element-(1) Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared)	Green ppt	
	(ii) Add few drops of dil. H ₂ SO ₄	Solution does not turn green or blue.	N-absent
4	Test for functional group FeCl ₃ test-Aqueous solution of the compound + neutral FeCl ₃ solution	Solution turns green	Phenolic (Ar-OH) group is present.
5	Liebermann Nitroso Test (i) Compound + solid NaNO ₂ + dilute H ₂ SO ₄ (ii) Dilute the solution with H ₂ O (iii) Above red solution + NaOH solution	[Blue coloured solution is obtained.] [Solution becomes red] [The colour of the solution again becomes blue.]	Phenolic group is present.

Result – Phenolic (Ar-OH) group is present in the given organic compound.

Structural Formula -



Experiment No. 8

Exercise – Identify the functional group present in the given organic compound.

Physical Properties –


1. Physical State – liquid
2. Colour – Colourless
3. Odour – Pungent smell
4. Solubility in water – Soluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound-Put a drop of liquid on moist blue and red litmus paper separately	No action	Organic compound is neutral.

S.N	Experiment	Observations	Inference
2	Burning Test – Put a crystal of the compound on glass rod and bring it in to the flame	The compound burns with non- sooty blue flame	Organic Compound is aliphatic.
3	Detection of element 1. Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared)	Green ppt	
	(ii) Add few drops of dil. H ₂ SO ₄	Solution does not turn green or blue.	N-absent
4	Test for functional group		
(i)	Compound + 2,4- dinitrophenyl hydrazine, shake	Yellow ppt is formed	Carbonyl (>C=O) group is present may be –CHO or > C=O group.
(ii)	Compound + sod. Nitroprusside+NaOH	No reaction	> C=O group is absent.
(iii)	Compound + Tollen's reagent warm in a water bath.	Silver mirror is formed.	–CHO group is present.
(iv)	Equal amount of Fehling's solution A and B in a test tube and add it to 1 ml given compound and heat.	Red ppt is formed.	–CHO present.
(v)	Schiffs reagent and the given organic liquid	Red colouration	Aldehyde group (–CHO) is present.

Result – Aldehyde group (–CHO) is present in the given organic compound.

Structural Formula -



Experiment No. 9 ✓

Exercise – Identify the functional group in the given organic compound.


Physical Properties –

1. Physical State – liquid
2. Colour – Colourless
3. Odour – Smell like nail polish remover
4. Solubility in water – Insoluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound: Put a drop of given liquid on moist blue and red litmus paper separately	No action	Organic compound is neutral.

2	Burning Test – Put a crystal of the compound on glass rod and bring it in to the flame	The compound burns with non- sooty blue flame	Organic Compound is aliphatic.
3	Detection of element-1. Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared)	Green ppt	
	(ii) Add few drops of dil. H ₂ SO ₄	Solution does not turn green or blue.	N-absent
4	Test for functional group		
(i)	Compound + 2,4- dinitrophenyl hydrazine, shake	Yellow ppt is formed	Carbonyl (>C=O) group is present may be –CHO or > C=O group.
(ii)	Totten's Test Organic compound + Tollen reagent	No reaction	–CHO group absent.
(iii)	Nitroprusside test – Organic compound + sodium nitroprusside solution + NaOH	A wine red colouration.	> C=O group is present.
(iv)	m-dinitrobenzene test – Organic compound + powdered m-dinitrobenzene + dil. NaOH in excess, shake well.	A violet colour appears	(>C=O) ketonic group is present.

Result – The given organic liquid has Ketonic (>C=O) functional group.



Experiment No. 10

Exercise – Identify the functional group in the given organic compound.


Physical Properties –

1. Physical State – liquid
2. Colour – Brownish red (Colourless when pure)
3. Odour – Pungent smell
4. Solubility in water – Partly soluble in water.

S.N	Experiment	Observations	Inference
1	Nature of the compound - Put a drop of given liquid on moist blue and red litmus paper separately	Red litmus paper turns blue	The compound is basic in nature.
2	Burning Test – Bring a glass rod dipped in the given liquid into the flame.	Compound burns with sooty flame	Given compound is aromatic
3	Detection of element Lassaigne solution + NaOH + Fe SO ₄ (Freshly prepared) add few drops of dil. H ₂ SO ₄	Green ppt A blue colour appears	N- is present.
4	Test for functional group		
(i)	Organic compound + few drops of conc. HCl + solid NaNO ₂ . Cool the solution in ice bath. Add alkaline B naphthol is cold.	Orange red dye is formed.	Aromatic amino (Ar-NH ₂) group is present.

Result – Amino (-NH₂) group is present attached to aromatic ring i.e. Aromatic amine (Ar-NH₂) is present.

Structural Formula -



Experiment No.11

Exercise – Analyze the given mixture for one acid and one basic radical by dry and wet tests.

S.N	Experiment	Observations	Inference
	Test for acid Radicals.		
1.	Mixture + dil. H_2SO_4	No reaction	Dilute add group is absent.
2	Mixture + conc. H_2SO_4 . Heat	Yellow brown fumes with pungent smell evolved	Conc. Acid group is present may be Br^-
3	Confirmatory test Silver Nitrate Test – Sodium Carbonate extract + dil HNO_3 + $AgNO_3$ solution.	Light yellow ppt is formed.	
4	Add NH_4OH in excess into above ppt.	Yellow ppt is partly dissolve.	Br^- confirmed.
5	CCl_4 layer test – Sodium carbonate extract + dil. HNO_3 + CCl_4 + Cl_2 water drop by drop and shake	CCl_4 layer turns yellowish orange.	Br^- confirmed.
	Test for basic Radical		
1	Mix + $NaOH$ solution. Heat	Smell of NH_3	Zero group is present may be NH_4^+
	Confirmatory Test		
2	Bring a glass rod dipped in HCl near the mouth of the test tube.	White fumes of NH_4Cl formed	NH_4^+ confirmed.
3	Bring a piece of filter paper moistened with Nessler's. Reagent near the mouth of the test tube.	Filter paper turns reddish brown	NH_4 confirmed
4	Bring a piece of filter paper moistened with $CuSO_4$ solution near the mouth of the test tube.	Filter paper turns deep blue	NH_4^+ confirmed

Result – The given mixture contains bromide (Br^-) acid radical and Ammonium (NH_4^+) as basic radical.

Experiment No. 12

Exercise : Identify one acid and one basic radical in the given mixture.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1.	<u>Tests for acid radicals</u> Mixture + dil H_2SO_4	Smell like vinegar is given out	Dilute acid group is present may be CH_3COO^-
2	Confirmatory Test Palm test – Take the mixture on palm, add it to solid $H_2C_2O_4$, $2H_2O$ and 2-3 Drops of water. Rub it and smell	Smell of vinegar	CH_3COO^- confirmed
3	Feric Chloride test-Sodium carbonate extract + neutral $FeCl_3$ solution	Solution turns deep red	CH_3COO^- confirmed
4	Ester test – Mixture + C_2H_5OH + Few drops and conc. H_2SO_4 and heat.	A pleasant fruity smell.	CH_3COO^- confirmed
Test for Basic Radical			
1	Mixture + $NaOH$ solution & Heat.	No smell of NH_3	Zero group absent.
2	Original solution + dil. HCl	White ppt	I group present may be Pb^{2+} , Ag^+ , Hg_2^{2+}
3	White ppt + water & boil	Ppt is soluble but reappear on cooling	May be Pb^{2+}
4	Divide the hot solution of the ppt into two parts.		
(i)	To part I add KI solution	Yellow part	Pb^{2+} confirmed
(ii)	To part II add K_2CrO_4 Solution Dry Test : Charcoal Cavity Test mixture + Na_2CO_3 +2-3 drops of water in charcoal cavity. Heat it in reducing flame.	Yellow part A soft shining metallic bead is formed. It leaves mark on paper.	Pb^{2+} confirmed Pb^{2+} confirmed

Result : - Acetate (CH_3COO^-) and lead (Pb^{2+}) are present in the given mixture as acid and basic radicals.

Experiment No. 13

Exercise : Analyze the given salt for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Salt + dil H_2SO_4 and heat	No reaction	Dilute acid group is absent.
2	Salt + conc. H_2SO_4 and heat	No reaction	Conc. Acid group is absent.
3	Sodium carbonate extract + dill $HCl+BaCl_2$ solution	White ppt is obtained	General groups is present may be SO_4^{2-}
4	Confirmatory test - Divide the ppt in two parts. (i) To part I add conc. HCl & Heat. (ii) To part II add conc. HNO_3 and heat.	Ppt remains insoluble. Ppt remains insoluble.	SO_4^{2-} confirmed SO_4^{2-} confirmed.
Test for Basic Radical			
1.	Salt + $NaOH$ solution & Heat.	No smell of NH_3	Zero group absent.
2	Original solution + dil. HCl	No ppt is obtained.	i group absent.
3	Pass H_2S gas into the above solution.	Black ppt is obtained	II group present may be Pb^{2+} or Cu^{2+}
4	Dissolve the ppt in 50% HNO_3	Black ppt is dissolves	May be Pb^{2+} or Cu^{2+}
5	Divide the solution in two parts. (i) Part I + C_2H_5OH + dil. H_2SO_4 (ii) Part II + NH_4OH in excess.	No ppt Deep the colouration	Pb^{2+} absent may be Cu^{2+} Cu^{2+} confirmed
6	Deep blue solution + CH_3COOH till blue colour disappear. Now add $K_4[Fe(CN)_6]$ solution.	Chocolate brown ppt is formed.	Cu^{2+} confirmed
7.	Dry test - Perform charcoal cavity test.	Red bead is formed.	Cu^{2+} confirmed

Result : - In the given salt sulphate (SO_4^{2-}) and copper (Cu^{2+}) radicals are present.

Experiment No. 14

Exercise : Analyze the given salt for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Mixture + dil H_2SO_4 and heat	No reaction	Dilute acid group is absent.
2	Salt + conc. H_2SO_4 and heat	No reaction	Conc. Acid group is absent.
3	Sodium carbonate extract + dill $HCl + BaCl_2$ solution	White ppt is obtained	General groups is present may be SO_4^{2-}
4	Confirmatory test – Divide the ppt in two parts. (i) To part I add conc. HCl & Heat. (ii) To part II add conc. HNO_3 and heat.	Ppt remains insoluble. Ppt remains insoluble.	SO_4^{2-} confirmed SO_4^{2-} confirmed.
	Test for Basic Radical		
1.	Mixture + $NaOH$ solution & Heat.	No smell of NH_3	Zero group absent.
2	Original solution + dil. HCl	No ppt	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2 Pb$ solution.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 + add solid NH_4Cl and then NH_4OH in excess after cooling.	Reddish brown ppt is obtained.	III group present may be Fe^{3+}
6	Confirmatory Test Dissolved the ppt in dil. HCl and divide it in two parts. (i) Part I + $K_4[Fe(CN)_6]$ solution. (ii) Part II + $KCNS$ solution Dry test borax bead test Perform borax bead test.	Deep blue colouration. Blood red colouration. Yellow bead.	Fe^{3+} confirmed. Fe^{3+} confirmed. Fe^{3+} confirmed.

Result : - The given mixture contains sulphate (SO_4^{2-}) and iron (Fe^{3+}) radicals.

Experiment No. 15

Exercise : Analyze the given mixture for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Mixture + dil H_2SO_4 and heat	No reaction	Dilute acid group is absent.
2	Mixture + conc. H_2SO_4 and heat	Brown coloured gas with pungent smell is given out.	Conc. Acid group is present. May be NO_3^-
3	Add some copper turning to the test tube.	Brown gas becomes dense.	Br- absent. May be NO_3^-
4	Confirmatory test – Neutralize Sodium carbonate extract with CH_3COOH . Add to it freshly prepared $FeSO_4$ solution and then conc. H_2SO_4 by the side of the test tube.	A brown ring is formed at the junction of two liquid due to the formation of $FeSO_4NO$	NO_3^- confirmed.
	Test for Basic Radical		
1.	Mixture + NaOH solution & Heat.	No smell of ammonia	Zero group absent.
2	Original solution + dil. HCl	No ppt.	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2 Pb$ solution.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	White gelatinous ppt is obtained.	III group present may be Al^{3+}
6	Confirmatory Test Dissolved the white ppt in dil. HCl and divide it in three parts. (i) Part I + solid NH_4Cl heat. (ii) Part II NaOH solution (iii) Lake test. Part III + 2 drops of blue litmus and then add NH_4OH drop by drop. Dry Test-Charcoal cavity cobalt nitrate test.	White ppt reappears. White ppt, which is soluble in excess if NaOH. Blue ppt floats on the solution. Blue mass is obtained.	Al^{3+} confirmed. Al^{3+} confirmed. Al^{3+} confirmed. Al^{3+} confirmed.

Result : - Nitrate (NO_3^-) and Aluminium (Al^{3+}) are present in the given mixture.

Experiment No. 16

Exercise : Identify one acid and one basic radical from the given mixture by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Mixture + dil H_2SO_4 + Heat	Colourless gas with smell of vinegar is evolved	Dilute acid group is present may be CH_3COO^-
2	Confirmatory test – Palm test : Mixture + solid oxalic acid on palm then add 2-3 drops of water. Rub it and smell.	Smell of vinegar	CH_3COO^-
3	$FeCl_3$ solution : Neutralize sodium carbonate extract with dil. HCl and add $FeCl_3$ solution.	Dark red colouration	CH_3COO^- confirmed
4	Ester Test : Mixture + C_2H_5OH + few drops of conc. H_2SO_4 . Heat	Pleasant fruity smell of ester formed	CH_3COO^- confirmed
	Test for Basic Radical		
1.	Mixture + NaOH solution & Heat.	No smell of ammonia	Zero group absent.
2	Original solution + dil HCl	No ppt.	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2 Pb$ solution.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	No ppt	III group absent.
6	Pass H_2S gas in the above solution.	A white ppt is obtained	IV group present may be Zn^{2+}
7	Dissolve the ppt in dil. HCl and divide the solution in two parts.		
(i)	Part I + $K_4[Fe(CN)_6]$ solution	Bluish white ppt is obtained	Zn^{2+} confirmed
(ii)	Part II + NaOH solution drop by drop	White ppt is formed which is soluble in excess of NaOH	Zn^{2+} confirmed
8	Dry test : Perform charcoal cavity test with cobalt nitrate.	Green residue is left.	Zn^{2+} confirmed

Result : - The given mixture contains acetate (CH_3COO^-) and Zinc (Zn^{2+}) in it.

Experiment No. 17

Exercise : Analyse the given salt for one acid and one basic radical by dry and wet tests

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Salt + dil H_2SO_4 and heat. Salt + conc. H_2SO_4 & heat. Bring a glass rod dipped in NH_4OH near the mouth of the test tube.	No reaction	Dilute acid group is absent.
2		Colourless pungent smelling gas evolved. White fumes formed.	Conc. Acid group is present. May be Cl^- Cl^- confirmed.
3	Confirmatory test – Chromyl Chloride Test Salt + solid $K_2Cr_2O_7$ + conc. H_2SO_4 in a dry test tube and heat. Pass the red coloured gas in $NaOH$ solution. To the yellow solution add CH_3COOH and $(CH_3COO)_2Pb$ solution.	Dark red pungent smelling gas is evolved. Solution becomes yellow. A yellow ppt is formed.	Cl^- confirmed. Cl^- confirmed.
4	Silver Nitrate Test : Sodium Carbonate extract + dil. HNO_3 + $AgNO_3$ solution. Add excess of NH_4OH in the ppt.	A white ppt is formed. The ppt is soluble.	Cl^- confirmed. Cl^- confirmed.
Test for Basic Radical			
1.	Salt + $NaOH$ solution & Heat.	No smell of ammonia.	Zero group absent.
2	Original solution + dil. HCl	No ppt.	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2Pb$.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	No ppt	III group absent.
6	Pass H_2S gas in the above solution.	No ppt	IV group absent.
7	Boil off H_2S from the above solution and concentrate it. Now add excess of NH_4OH and saturated solution of $(NH_4)_2CO_3$.	A white ppt is obtained.	V group present. May be Ba^{2+} , Sr^{2+} and Ca^{2+}
8	Dissolve the white ppt in hot CH_3COOH and divide it in to three parts.		
(i)	To part I add $K_2C_2O_4$ solution. Dry test – Perform flame test.	A yellow ppt is obtained Apple green flame.	Ba^{2+} confirmed. Ba^{2+} confirmed.

Result - The given mixture contains chloride Cl^- and Barium Ba^{2+} radical.

Experiment No. 18

Exercise : Analyse the given salt for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Salt + dil H_2SO_4 and heat.	No reaction	Dilute acid group is absent.
2	Salt + conc. H_2SO_4 & heat. Bring a glass rod dipped in NH_4OH near the mouth of the test tube.	Colourless pungent smelling gas evolved. White fumes formed.	Conc. Acid group is present. May be Cl^- Cl^- is confirmed.
3	Confirmatory test – Chromyl Chloride Test Salt + solid $K_2Cr_2O_7$ + conc. H_2SO_4 in a dry test tube and heat. Pass the red coloured gas in NaOH solution. To the yellow solution add CH_3COOH and $(CH_3COO)_2Pb$ solution.	Dark red pungent smelling gas is evolved. Solution becomes yellow. A yellow ppt is formed.	Cl^- confirmed. Cl^- confirmed
4	Silver Nitrate Test : Sodium Carbonate extract + dil. HNO_3 + $AgNO_3$ solution. Add excess of NH_4OH in the ppt.	A white ppt is formed The ppt is soluble.	Cl^- confirmed Cl^- confirmed.
Test for Basic Radical			
1.	Salt + NaOH solution & Heat.	No smell of ammonia.	Zero group absent.
2	Original solution + dil. HCl	No ppt.	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2Pb$.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	No ppt	III group absent.
6	Pass H_2S gas in the above solution. Boil off H_2S from the above solution and concentrate it. Now add excess of NH_4OH and saturated solution of $(NH_4)_2CO_3$. Dissolve the white ppt in hot CH_3COOH and divide it in to three parts.	No ppt A white ppt is obtained.	IV group absent. V group present may be Ba^{2+} , Sr^{2+} or Ca^{2+} .
7			
8			
(i)	To part I add K_2CrO_4 solution.	No ppt	Ba^{2+} absent.
(ii)	To Part II add $(NH_4)_2SO_4$ solution	A white ppt is obtained.	Sr^{2+} confirmed.
(ii)	Dry test – Perform flame test with the salt or the ppt.	Consistent crimson red flame	Sr^{2+} confirmed.

Result : - Chloride (Cl^-) and strontium (Sr^{2+}) radicals are present in the given salt.

Experiment No. 19

Exercise : Analyze the given mixture for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Test for acid radical Salt + dil H_2SO_4 and heat.	Brisk effervesce with evolution of colourless and pungent smell like burning sulphur gas.	Dilute acid group is present may be SO_3^{2-}
2	Pass the above gas in lime water.	Lime water turn milky.	SO_3^{2-} confirmed.
3	Pass more of the gas in milky solution.	Milkyness disappears.	SO_3^{2-} confirmed.
4	Bring acidity $K_2Cr_2O_7$ clipped filter paper on the mouth of the test tube.	Filter paper turn into green.	SO_3^{2-} confirmed.
Test for Basic Radical			
1.	Salt + NaOH solution & Heat.	No smell of NH_3	Zero group absent.
2	Original solution + dil. HCl	No ppt.	I group absent.
3	Pass H_2S gas into the above solution.	No ppt.	II group absent.
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2 Pb$.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	No ppt	III group absent.
6	Pass H_2S gas in the above solution.	No ppt	IV group absent.
7	Boil off H_2S from the above solution and concentrate it. Now add excess of NH_4OH and saturated solution of $(NH_4)_2CO_3$.	A white ppt is obtained.	V group present. May be Ba^{2+} , Sr^{2+} & Ca^{2+}
8	Dissolve the white ppt in hot CH_3COOH and divide it in to three parts.		
(i)	To part I + K_2CrO_4 solution.	No ppt.	Ba^{2+} absent.
(ii)	To Part II + $(NH_4)_2SO_4$ solution.	No ppt.	Sr^{2+} absent.
(iii)	To Part III + $(NH_4)_2C_2O_4$ solution.	White ppt.	Ca^{2+} confirmed.
9.	Dry Test – Perform flame test.	Brick red flame.	Ca^{2+} confirmed.

Result : - The given mixture contains _____ and calcium (Ca^{2+}) radicals.

Sulphite (SO_3^{2-})

Experiment No. 20

Exercise : Analyse the given salt for one acid and one basic radical by dry and wet tests.

S.N	EXPERIMENT	OBSERVATION	INFERENCE
1	Tests for acid radicals Salt + dil H_2SO_4 and heat.	Brick effervesence with the evolution of colour less and odourless gas.	Dilute acid group is present. May be CO_3^{2-}
2	Confirmatory test – Pass the evolved gas in lime water.	Lime water turns milky.	CO_3^{2-} Confirmed.
3	Pass more of the gas in milky solution.	Milkyness disappears.	CO_3^{2-} Confirmed.
	Test for Basic Radical		
1.	Salt + NaOH solution & Heat.	No smell of NH_3	Zero group absent.
2	Original solution + dil. HCl	No ppt .	I group absent.
3	Pass H_2S gas into the above solution.	No ppt	II group absent
4	Boil the above solution to remove H_2S gas completely. Test it by a piece of filter paper moistened with $(CH_3COO)_2 Pb$.	Filter paper does not turn black.	H_2S is removed completely.
5	Now boil the above solution with 2 drops of conc. HNO_3 add solid NH_4Cl and then NH_4OH in excess after cooling.	No ppt	III group absent.
6	Pass H_2S gas in the above solution.	No ppt	IV group is absent.
7	Boil off H_2S from the above solution and concentrate it now add excess of NH_4OH and $(NH_4)_2 CO_3$	No ppt	V group absent.
8	To the solution obtained above add disodium hydrogen phosphate. (Na_2HPO_4)	White ppt is formed. Which is soluble in dil. HCl	VI group present may be Mg^{2+}
9	Solution of white ppt in dil. HCl + NaOH + Titan yellow. Dry test – Perform charcoal cavity test with Cobalt Nitrate.	Pinkish red ppt is obtained. Pink residue is left in the cavity.	Mg^{2+} confirmed. Mg^{2+} confirmed.

Result : - The given salt contains carbonate (CO_3^{2-}) and Magnesium (Mg^{2+}) radicals in it.